Not for Use in USA

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Not for Use in USA

Application and safety instructions for the MEISINGER

EN

Micro Screw System Micro Screw System Basic

Screw System for safe Stabilization of Cortical Bone Graft Transplants ‘Olsberger Konzept’ by Prof. Dr. Fouad Khoury

Art.-No.: BMS00, BMSBA

The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1,0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

Indications for Use

The Bone Management® Micro Screw Systems are indicated for the treatment of bone fractures, especially for the fixation of transplanted bone blocks during the augmentation process in the oral cavity and maxillomandibular surgical field. Note: Micro Screws are not intended to remain in the body permanently. After they have fulfilled their supportive function such as is the case after healing of a transplant, or healing of a fracture, for example, they need to be removed completely.

Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Placement of the pilot hole through bone cylinder and recipient bone


For the application of the Micro Screws, a pilot hole through bone graft and recipient bone is necessary. This is done using a twist drill with a 0.2 mm smaller diameter:

- 1,0 mm Micro Screw ≙ twist drill MSPB1 or initial bur 202RF 008 (Ø 0.8 mm)
- 1.2 mm Micro Screw ≙ twist drill MSPB2 or initial bur 202RF 010 (Ø 1.0 mm)

The respective drilling should take place intermittently and under constant cooling with sterile, physiological saline solution. Since the Micro Screws have a self-tapping thread due to specially shaped thread flanks, no expansion drilling (gliding hole) to the outer diameter of the screws is necessary.

Selection of the screwdriver and picking up of the Micro Screws

Depending on the selected procedure (machine or manual screw insertion), the appropriate screwdriver must be selected. The screwdrivers MSSDE (screwdriver basic) or MSSDM (optionally with claw MSGSD) are used for manual insertion. For machine insertion with the angle piece, the MSSDW, screwdriver RA, is used.

Fixation of the bone graft

Using the selected screwdriver, a Micro Screw of appropriate length is picked up and the bone graft is attached to the recipient bone. The Micro Screw can be inserted into the bone without thread cutting, as it has a self-tapping thread design. The Micro Screw is tightened in a controlled manner until the bone cylinder is firmly fixed in the site. The bone must not be over stressed. In addition, during the insertion process, the axis specified by the pilot hole must be observed to prevent the screw from tilting. If a safe fixation of the graft with one screw is not possible, fixation with at least 2 screws is recommended to prevent any rotation.

Removing the Micro Screws

Taking into account the patient-specific healing time and after clinically and radiographically confirmed complete bone healing as well as functional restoration, the Micro Screws are to be removed counter-clockwise using a screwdriver from the Micro Screw System. The removal of the Micro Screws is also necessary in the case of screw loosening, infection or weakening of the bone tissue due to the inserted bone screw.

Sterilisation of the Micro Screws

The instruments and Micro Screws are non-sterile when delivered and must be processed (cleaned / disinfected / sterilised) according to validated instructions for processing before use. For proper processing of the articles, please observe the instructions for processing (cleaning, disinfecting and sterilising) of medical devices produced by Hager & Meisinger GmbH.

Instructions for disassembly and assembly of the screwdriver with claw

Screwdriver ①, Claw ②, Groove 1 ③, Groove 2 ④

1. Attach the claw ② to the screwdriver blade ①, as shown. Push the claw ② until it snaps into groove 1 ③, open position).

2. Now the screwdriver with claw is ready to use. To secure a screw, push the claw ② in the direction of the screwdriver blade over the screw until it snaps into groove 2 ④, closed position). For Disassembly push and remove the claw ② into the direction of the screwdriver blade.

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Disposal
When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. Dual Waste System).

Contraindications
Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

Absolute contraindications

- Active infections in or near the area to be augmented as well as local or systemic pathological processes (e.g. symptoms such as fever, local inflammation, abscesses)
- Insufficient bone volume (quality / quantity)
- Diseases or medication that impair bone metabolism

Relative contraindications

- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Known allergies and/or foreign body reaction to the alloy composition of the Micro Screws
- Systemic and/or metabolic diseases or medical treatments that lead to progressive deterioration of the bone (e.g. cortisone, immunosuppressants, bisphosphonates, diabetes, etc.)
- Drug and alcohol abuse
- Lack of patient compliance
- Poor circulation
- Hard physical labor or active sports
- Mental condition, which can lead to disregard of the medical order
- Highly atrophic jaw

IMPORTANT: Attention must be paid to the protection of the anatomical structures (safety distance min. 2 mm) as well as to the course of the adjacent teeth / tooth roots (risk of damage, infection, resorption). To avoid infection and resorption, care should be taken to avoid contact between the graft and adjacent teeth. Mechanical stability is essential for successful bone healing.

SEPARATE NOTES FOR SCREWS: The bone screws are compatible with the instruments of the Micro Screw System and Micro Screw System Basic. The screws may only be used with these system instruments. When selecting the screw length as well as the screw placement, the evaluation of the available bone by means of X-ray diagnostics must be taken into account. The screw head should be countersunk as flat as possible in the bone block to prevent soft tissue perforation or dehiscence. Excessive physical activity as well as trauma impairing the augmentation site can lead to premature failure of the screws due to loosening or breaking. The screws may only be used in accordance with the application and safety instructions. Any kind of modification in the procedure or on the product can lead to impairment of the success of the treatment or damage to the product. The screws are a disposable product and must therefore not be reprocessed and reused after use. Even with screws that have already been used but are visually undamaged, the previous mechanical stress can cause microscopic damage, which can lead to fatigue and product failure. The Micro Screws have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Micro Screws in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

RISK OF ASPIRATION: The user must ensure that the screws are secured against aspiration both during the attachment and during its removal.

Any serious-incident that has occurred in relation to the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

General instructions
Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Manufacturer	Sales partner	Item number	Medical device	Batch number	Do not re-use	Non-sterile	Follow the instruction for use	Under US Federal law, this device may only be sold to or on behalf of trained medical professionals.	Caution

CE / CE 0044
Conformity mark

Do not use if the packaging is damaged

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