Meisinger

BONE MANAGEMENT® CATALOGUE

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When developing our surgical systems and instruments, we are not only concerned with the highest quality of "Made in Germany", but also with the correct application of these systems and instruments. Since the professional preparation of the bony implant site is essential for successful implant placement and management, MEISINGER has developed a multi-level training concept in cooperation with leading scientists and practitioners. With the Bone Management® Academy, both experienced users and beginners have the opportunity to deal intensively with the topics of bone augmentation, implantology, content related to practice management as well as providing the chance to further their education in this field.

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
Our main motivation is on the continuously increasing added value for participants.

The focus is therefore not on the systems developed by MEISINGER, but on professional, practical, and patient-oriented techniques that can be learned in the courses. High-quality and precisely manufactured products play an essential role in positive and successful patient outcomes. However, the correct and proficient application of these tools is the famous "key to success" and thus also the best way to a satisfied patient.

That is why we do not shy away from our competition, but deliberately involve other dental companies in order to enable our participants to form a balanced perspective, and to discuss and compare competing treatment methods. If the participants find certain techniques to be good, the products can firstly be intensively tested and practiced in the workshops. Additional work shadowing and supervision offer further opportunities to consolidate what has been learned and thus make better informed decisions.

The Young dental experts camps, which take place several times a year, also offer orientation aid for students, interns, and residents in their future specialization choice. Aspiring surgeons can also acquire solid basic knowledge in the various courses and workshops. The right suturing and cutting techniques are just as important as the right understanding of soft and hard tissue. Various augmentation techniques can also be learned and practiced within the curriculum. The Bone Management® Academy also includes complex techniques such as the Carrot and Shell techniques developed by Prof. Dr. Fouad Khoury or the Sausage technique according to Prof. Dr. Istvan Urban.

SITES / LOCATIONS

HAGER & MEISINGER GMBH :

-  MEISINGER GERMANY
-  MEISINGER USA
-  MEISINGER FRANCE SARL
-  MEISINGER ASIA CO. LTD.

SITES / LOCATIONS

MYPLANT GMBH :

-  MYPLANT UK & IRELAND LTD.
-  MYPLANT TÜRKIYE AS
-  MYPLANT ITALIA S.R.L.

SOMMAIRE/CONTENT

08 EXTRACTION DE RACINE

08-09 CBE00 | Benex®-Control

10 AUGMENTATION

10-11 CCR00 | Crest-Control
12-13 CSP12 | Split-Control 12 mm
14-15 CSPPL | Split-Control Plus
16-17 BSPPR | Split-Control Professional

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64 BSK03 | Surgical Kit 3
64 BSK04 | Surgical Kit 4
65 ALV18 | Alveoplasty Surgical Kit
66 IPK02 | Implant Preparation Kit with Stop
67 BDS00 | Drill-Stop-Control (BL+TL)
68 BDS02 | Drill-Stop-Control (BLT)
69 BDSMP | Drill-Stop-Control (myplant two)
70-71 BGS00 | Guided-Drill-Stop-Control
72 BNA00 | 3D-Navigation-Control
73 BNAPR | 3D-Navigation-Control Professional

76 PERIIMPLANTITES

74 2575 | Periimplantitis Kit
75 DEG00 | Degranulation Kit

18 BONE MANAGEMENT® MASTER-LINE

20 BMCBA | Master-Core Basic
21 BMCPR | Master-Core Professional
22-23 BKM00 | Master-Mill
24 BMPBA | Master-Pin-Control Basic
25 BMP00 | Master-Pin-Control
26-27 Master-Pin-Control Plus

28 BONE MANAGEMENT® KHOURY-LINE

30-31 BTE00 | Trephine Ejection Kit
32 BMSBA | Micro Screw System Basic
33 BMS00 | Micro Screw System
34-35 BMS00 | Micro Screw System Plus
36-37 BOST0 | Ost-Tray
38 BSINO | Sinus-Tray
39 BSOTO | Soft-Tissue-Tray
40-41 BS003,BS004 | SAFESCRAPER® Twist

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42-43 CTRBA | Transfer-Control Basic
44-45 CTR00 | Transfer-Control
46 CTRIO | Transfer-Ring-Control I
47 BBR00 | Transfer-Ring-Control II
48-49 CTRPL | Transfer-Control Plus
50 BTX00 | Screw System TX
51 BTXPR | Screw System TX Professional
52 7120 | Threphine Basic Kit
53 7121 | Threphine Kit
54 7122 | Threphine Kit Long
55 7140 | Punch Basic Kit- Saw basic kit

56 SINUSLIFT

56 CCLBA | Crestal-Lift-Control Basic
57 BCL00 | Crestal-Lift-Control
58 BTL00 | External-Lift-Control
59 CSL00 | Internal-Lift-Control
60-61 BLIPL | Lift-Control Plus
62 BSKSL | Surgical Kit 1
63 BSK02 | Surgical Kit 2

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Benex®-Control is outstandingly suitable for the controlled and safe extraction of dental roots and dental root fragments, so the conventional extraction method using forceps and leverage can be avoided. Thanks to the innovative construction of the extractor, the root can be easily removed in a particularly controlled manner without damaging bone or soft tissue. To guarantee an optimal vertical transmission of force, there must be a flat area providing a supporting surface to the extractor. The quadrant support was developed specifically to provide this optimal support while also protecting existing dentition and providing support for the extractor in case of edentulous circumstances.

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
* Plaque de support Teflon® inclus incl. straight Teflon® bite plate

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AT A GLANCE

- Safe and controlled extraction of roots
- Possibility of immediate implant placement
- Maximum preservation of the vestibular bone wall due to vertical extraction
- Optimal protection of bone and soft tissues
- Applicable to small root fragments
- Osteotomy can be prevented
- High patient acceptance

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

Crest-Control is a Horizontal Bone Splitting System which allows for a predetermined and minimally invasive widening of the horizontally resorbed alveolar ridge particularly in the distal lower jaw. With the well-proven Horizontal Spreaders, the alveolar ridge can be expanded up to a maximum of 5 mm, allowing an implant of every common system to be inserted in the widened segment.

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AT A GLANCE

- Controlled widening of horizontally resorbed bone
- Effective in increasing the width of narrow alveolar ridges
- Preparation of the bone for all common implant systems
- Minimally invasive treatment concept
- Gentle widening of the alveolar ridge

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
Art.-No. CSP12

Split-Control 12 mm is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

Split-Control 12 mm provides an extensive range of instruments for a controlled dilatation with spreaders with a length of 12 mm.

	Pointeur Initial bar	Scies à revêtement diamanté Diamond coated saws	Instruments d'ostéotomie Osteotomy instruments	Foret Pilote Pilot burs	Foret d'expansion 12 mm Expansion burs 12 mm	Clé à cliquet Ratchet				
Ref.	186RF	231DC*	231DC*	HM33IL**	HM254E**	B1001	E1001	A2003	D2003	CARA4
Tige ¹	204	204	204	205	205	206	206	205	205	-
Taille ²	018	070	130	010	012	013	020	023	030	-
Longueur mm	12.0	0.3	0.3	5.5	6.0	18.0	18.0	12.0	12.0	84.0
Δ	-	-	-	-	-	-	-	1.30	2.03	-
□	1.8	7.0	13.0	1.0	1.2	1.3	2.0	2.3	3.0	-

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AT A GLANCE

- Controlled spreading of horizontally resorbed bone
- Gentle bone condensing with the aid of non-cutting Archimedes screws
- Increased primary stability of the inserted implants
- Effective in increasing the width of narrow alveolar ridges
- Preparation of the bone for all common implant systems

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
Split-Control Plus is a combination of the proved Split-Control and Crest-Control systems, so it offers an extensive range of instruments for controlled bone spreading and condensing. It contains both specially designed screw-like condensing and spreading instruments (Spreader) and horizontal spreaders that allow for controlled and standardized spreading of horizontally resorbed bone. Due to the special geometry of the screw-like spreaders, cancellous bone is gently condensed while the Horizontal Spreaders allow for a dilatation up to a width of 5 mm, so the bone is optimally prepared for the insertion of every common implant.

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

AT A GLANCE

- Combined System of Crest- and Split-Control
- Possibility of combining bone splitting and spreading techniques
- Horizontal spreaders for gentle widening of the alveolar ridge by up to 5 mm
- Non-cutting Archimedes screws for gentle bone condensing
- Preparation of the bone for insertion of all common implant systems

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Split-Control Professional is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased. To be able to respond more precisely in specific indications, the Split-Control Professional contains a broad range of spreader sizes.

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
BONE MANAGEMENT® COMPETENCE CENTER

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The Master-Core Basic is used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the system contains trephines with different lengths and diameters. Due to the different lengths, the trephines have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trephines' working parts enable additional depth control during the extraction. The black coating of the trephines allows for glare-free work and the optimal visibility of the depth markings.

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
AT A GLANCE

- Simple and safe extraction of bone cylinders
- Optimal protection of anatomical structures thanks to depth-limited trephines
- Depth markings on the trephines for additional depth control
- Extraction of bone cylinders of different diameters and lengths
- Black coating on the trephines for glare-free work


Ref.	DV229	DV229	BV025
Tige ¹	204	204	205
Taille ²	050	070	-
Longueur mm	2.0	2.0	25.5
	4.0	5.0	-
	4.7	5.7	-

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
The Master-Core Professional is used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the system contains trephines with different lengths and diameters. Due to the different lengths, the trephines have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trephines' working parts enable additional depth control during the extraction. The black coating of the trephines allows for glare-free work and the optimal visibility of the depth markings. To be able to respond more precisely in specific indications, the Master-Core Professional contains a broad range of instruments.

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

The Master-Mill is a surgical bone mill used for crushing autologous bone blocks and thus preserving small autologous bone chips. With the help of these bone chips, bone defects of all shapes can be optimally filled. The individual components of the Master-Mill are made of high quality, surgical stainless steel. All parts feature an optimal fit and a stable geometry, allowing the Master-Mill to be easy to handle with ease of assembly and intuitive use.

Master-
Line
par le
Pr. Istvan Urban

Art.-No. BKM00

Master-Mill (broyeur inclus)
Master-Mill (incl. Grinder)

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AT A GLANCE

- Crushing of autologous bone blocks
- Augmentation of bone defects of all shapes
- Individual components made of high quality, surgical stainless steel
- Stable geometry
- Easy and intuitive handling

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
Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The included pins consist of titanium alloy Grade 5 and have a groove which increases the surface of the pins. Due to the groove and a unique head-design, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.

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
AT A GLANCE

- No bending, no breakage - Safe Fixation of Membranes with Pins made of Titanium Alloy Grade 5
- Sharp tip and high stability of the pins allow for precise placement of the pins also in cortical bone and at an angle
- Easy Removal from the tray and placement of the pins due to the specially developed Pin-Holder
- Instruments for the placement of Bleeding Points for optimal healing conditions

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Ref.	186RF	203RF	203RF	MP10*
Tige ¹	204	205	205	-
Taille ²	018	006	008	-
Longueur mm	12.0	7.0	7.0	3.65
Δ	-	-	-	0.87
+ □	1.8	0.6	0.8	0.95
Ø tête mm	-	-	-	2.51

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Twist Drill with Stop for decorticating

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
Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The Master-Pin-Control Plus contains two different types of pins. While one pin has a particularly flat head design and is therefore, for example, ideal for use under thin soft tissue, the other pin has a larger head design that makes it particularly easy to find the pins under the soft tissue after successful healing. All included pins consist of titanium alloy Grade 5 and have a groove which increases the surface of the pins. Due to the groove, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.


For the simple and safe placement of bleeding points in the cortical bone, the Master-Pin-Control Plus system also contains the twist drill 203S, which has an automatic depth stop.

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AT A GLANCE

- Will not bend or break - Secure attachment of membranes with grade 5 titanium pins.
- 2 different types of pins for optimal adaptation to all types of soft tissue.
- These particularly sharp and stable pins allow precise pin placement, even in dense cortical bone, and even at an angle.
- Easy removal of the pins from their holders thanks to the specially developed gripper.
- 203S stop drill for safe and easy creation of blood flow and for creating the best healing conditions.

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


The Trephine Ejection Kit was developed together with Prof. Dr. Fouad Khoury. It contains several unique ideas which turn this kit into something special: Two-part, internally cooled trephines allow for a safe and minimal invasive extraction of local bone cylinders. The inner cooling helps to prevent a damage of the bone tissue due to overheating, while the removable trephine working parts allow for an especially easy sampling of the cut out bone cylinders. Either with the aid of the ejection needle or the ejection instruments, which perfectly match with the diameter of the trephine working parts, the bone cylinders can be pushed out of the trephines. Furthermore, the included pre-drilling trephines allow for a safe application of the trephines. Due to the two-part trephines, the variety of different diameters and the optimal matching of the instruments, the Trephine Ejection Kit offers highest efficiency, flexibility and precision to the user for the extraction of bone cylinders.

Art.-No. BTE00

Bone core are harvested from the implant site minimal invasive with the pre- and main trephine in different diameter.

Bone core augmentation stabilized with micro-screw

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
AT A GLANCE

- Two-part trephines for a safe and minimal invasive extraction of bone cylinders
- Protection against overheating due to internally cooled trephines
- Innovative application
- Standardized diameters and precise laser markings of the trephines allow for a predictable extraction of defined bone cylinders
- Perfectly matched instruments care for highest precision and flexibility

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The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1.0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

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
The Micro Screw System Plus is a combined system of the Trephine Ejection Kit, which enables the minimally invasive extraction of local bone cylinders, and the Micro Screw System, which is used for the safe fixation and stabilization of bone grafts. The Micro Screw System Plus thus enables the complete workflow of the Bone-Core-Technique to be carried out. For the extraction of local bone cylinders, the system contains two-part, internally cooled trephines, which, thanks to the removable trephine working parts, enable the bone cylinders to be pushed out particularly easily with the help of special ejection instruments. For the subsequent fixation of the bone cylinders, the Micro Screw System Plus offers osteosynthesis screws with a diameter of 1.0 mm or 1.2 mm, as well as all the necessary instruments for setting the screws. The osteosynthesis screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

The combination of the two techniques in one system and the optimally matched instruments allow for the highest precision from the extraction to the fixation of autologous bone cylinders.

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
AT A GLANCE


- System for carrying out the complete workflow of the Bone-Core-Technique
- Optimally matched instruments for highest precision from extraction to fixation of autologous bone cylinders
- Two-part trephines for the safe and minimally invasive extraction of bone cylinders
- Safe fixation and stabilization of cortical bone grafts
- Osteosynthesis screws with diameters of 1.0 and 1.2 mm made from surgical stainless steel for high stability
- Self-tapping thread flanks of the screws allow the osteosynthesis screws to be screwed in without creating a gliding hole

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The Bone Management® Ost-Tray offers a perfect combination of specially configured manual surgical instruments. In addition to a well-thought-out selection of manual instruments, such as high-quality periodontal probes, curettes, scalpel blade holders, mouth mirrors, graduated mouth wedges, a curved irrigation cannula with bulb end and various scissors, needle holders, tweezers, clamps, forceps, retractors and tongue depressors, the kit contains an individually designed, strong Khoury-probe as well as raspatories and sharp spoons according to Prof. Khoury.

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

The new Bone Management® Sinus-Tray is ideal for performing the external sinus lift. The manual instruments designed with Prof. Khoury support the user in gently lifting the Schneiderian membrane and inserting the bone substitute material.

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
Inclusive: Bur Block BSCM0 (for Handpiece, 10 slots)

Instruments manuels
Manual instruments

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
The new Bone Management® Soft-Tissue-Tray was developed with Prof. Khoury and is perfectly applicable in the context of soft tissue management. A special feature is the scalpel blade holder with swivel head included in the set, with which the scalpel angle can be freely adjusted depending on the processing area. This simplifies handling, especially in areas that are difficult to access.

Art.-No. BSOT0

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Instruments manuels
Manual instruments

Ref.	M1042	M1049	M1060	M1077	M1078	M1079	M1089	M1067
Nom	Porte-aiguille "Crile Wood", RH Needle Holder "Crile Wood", RH	Précelle "Cooley" Tweezers Cooley	Instruments de tunnelisation Tunneling Instrument	Porte-lame Scalpel Handle	Porte-lame Scalpel Handle	Porte-lame à tête pivotante Scalpel Handle with Swivel Head	Ciseau de chirurgie "Kelly" Surgical Scissors "Kelly"	Décolleur pour papille Papilla Elevator
Longueur mm	150.0	175.0	176.0	173.0	135.0	149.0	160.0	170.0

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With the Safescraper® Twist, autologous cortical bone chips can be obtained easily and safely. These are safely collected in the transparent chamber while the Safescraper® is in use - for optimal control of the amount of bone harvested. The bone material obtained can be used immediately thanks to the immediate mixture with blood and is suitable for filling up defects of all kinds.

Safescraper® TWIST Droit
Safescraper® TWIST
straight Version

Ref.

BS003

Safescraper® TWIST Angulé
Safescraper® TWIST Curve
angled Version

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
AT A GLANCE

- Ready to use
- Excellent cutting performance for a minimally invasive and patient-friendly procedure
- Straight and angled design for every situation
- Half-round blade with 160° cutting area for maximum efficiency on all bone surfaces
- Capacity of the chambers of 3.0 cm³ (straight version) or 2.5 cm³ (angled version)
- Safety mechanism prevents the chambers from opening during bone harvesting

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


Art.-No. CTRBA

Transfer-Control Basic permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in three diameters. They are perfectly adapted to each other so that the outer diameter of the abrasive bur and the size of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control Basic is an introductory system with reduced instrument set. It is suitable for the extraction of bone cylinders with diameters of 4.0 mm, 5.0 mm and 6.0 mm.

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AT A GLANCE

- Easy and safe transplantation of bone cylinders
- Intelligently coordinated instruments
- Perfectly fitting bone cylinders
- Precise work due to depth marks
- Optimally suitable for vertical and horizontal augmentations
- Creation of refreshed contact surfaces for rapid vitalization and wound healing

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
Transfer-Control permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in five diameters. They are perfectly adapted to each other so that the outer diameter of the ablative bur and the size of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control is suitable for the extraction of bone cylinders with diameters of 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm and 8.0 mm.

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


The most important condition for a safe healing of a bone graft is a congruent and fresh contact surface of the implant area. The Transfer-Ring-Control I System allows for this precondition for the healing of vertical bone grafts in a simple and controlled manner. This system offers a variety of various instrument sizes and provides high flexibility for the extraction of bone rings. Depending on requirements, various sizes of Initial Bur Trephines, Trephines and Ablative Burs are available. With the aid of the additional tungsten carbide bur, sharp bone edges can be smoothed precisely.

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

Transfer-Ring-Control II is a system for the bone ring technique, developed by Dr. Giesenhagen. It allows for one-stage proceeding of augmentation and implantation in case of resorbed jaw bone, if implant insertion alone is not possible because of inadequate spatial conditions. Due to the perfectly matched instruments, there is created a congruent bone and implant site, which allows for the treatment of three-dimensional defects with the help of prefabricated, allogenic bone rings. Because of the optimal fit between jaw bone, bone ring and inserted implant and the possibility of the one-stage proceeding, high primary stability is accomplished and healing time is significantly reduced.

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
Transfer-Control Plus permits precise and standardized transplantation of bone cylinders and rings for horizontal and vertical bone augmentation. Different types of bone burs allow for high flexibility for the extraction of bone grafts. The instruments are perfectly adapted to each other often producing a press fit between the cylindrical bone or ring and the prepared recipient site. This press fit can then be strengthened with a fixation screw which is sold separately (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound-healing. So a transplantable bone site is reached after only 3-4 months.

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


AT A GLANCE

- Combined system for easy and safe transplantation of bone cylinders and rings
- Intelligently coordinated instruments
- Optimally suitable for vertical and horizontal augmentations
- Precise work due to depth marks
- Creation of refreshed contact surfaces for rapid vitalization and wound healing

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

The Osteosynthesis System Screw System TX enables the universal use of fixation for autologous bone cylinders, bone blocks and bone shells for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

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
The Osteosynthesis System Screw System TX Professional enables the universal use of fixation for autologous bone cylinders, bone blocks and bone shells for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

The Screw System TX Professional contains additional screws in the lengths of 13.0 mm and 16.0 mm.

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The Trephine Basic Kit provides a selected set of high-quality trephine drills in different sizes. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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AT A GLANCE

- Selection of high quality trephines for the extraction of bone cylinders
- Various trephine sizes for every application
- Precise depth markings for optimal control
- High cutting performance for good handling

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The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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AT A GLANCE

- Broad selection of high quality trephines for the extraction of bone cylinders
- Various trephine sizes for every application
- Precise depth markings for optimal control
- High cutting performance for good handling

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The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 14 mm* and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

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
The Punch Basic Kit contains Biopsy Punches in five different diameters for the performance of tissue removal. They serve for the application in immediate implant placement, for taking mucosa grafts in the scope of soft tissue management or biopsy as well as for uncovering the gingiva after implantation. The Biopsy Punches create clean cut surfaces and care for minimal tissue traumatization.

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The Saw Basic Kit provides a selection of diamond saws in the three most important diameters. They are used for the application in bone spreading or splitting operations. Due to the various sizes, the set offers the clinician the possibility of having a proper selection of these saws during surgery.

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

Crestal-Lift-Control Basic is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcristal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control Basic is an introductory system with reduced instrument set. **CAUTION**: To prevent injury of any anatomical structure, the apical extra lengths of the CL001 of 0.58 mm must be considered.

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
Crestal-Lift-Control is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcristal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control provides an extensive range of instruments for the performance of a safe internal sinus lift.

CAUTION: To prevent injury of any anatomical structure, the apical extra lengths of the CL001 of 0.58 mm must be considered.

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
External-Lift-Control is a system for the easy performance of an external sinus lift. After lateral access through the facial wall of the maxillary sinus, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. Then, a suitable augmentation material can be inserted under the membrane. Due to the special CARBOCER[®]-coating, the hand-held instruments provide a reflection-free working environment.

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

Internal-Lift-Control is a system for the performance of an internal sinus lift. With the aid of the elevators, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 5 mm.

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The combined Lift-Control Plus system provides perfectly matched instruments for the performance of internal as well as external sinus floor elevation. For the internal sinus lift, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 5 mm. After lateral access through the facial wall of the maxillary sinus for the performance of the external sinus lift, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. After this, a suitable augmentation material can be inserted under the membrane.

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

The Surgical kit 1 contains instruments with handpiece shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. The large diameter of the instruments protects the membrane from unintentional perforation.

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AT A GLANCE

- Compilation of handpiece instruments for the lateral approach during the external sinus lift
- Round instruments protect the Schneiderian membrane from unintentional perforation
- Ideal for gradual, safe thinning of the facial wall of the maxillary sinus
- Round diamonds and carbide burs in three different diameters

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The Surgical Kit 2 contains instruments with right-angled shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. Especially the large diameter of the instruments protects the membrane from unintentional perforation.

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AT A GLANCE

- Compilation of contra-angle instruments for the lateral approach during the external sinus lift
- Round instruments protect the Schneiderian membrane from unintentional perforation
- Ideal for gradual, safe thinning of the facial bone wall
- Round diamonds and carbide burs in three different diameters

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The Surgical Kit 3 contains conical Tungsten Carbide Burs with friction grip shank of various sizes and tothing which allow for a broad field of application. The particularly easy cutting and efficient cutters have an extra-long shank for an optimal view on the working area and sharp tips for the safe application in bone.

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The Surgical Kit 4 contains a selection of conical Tungsten Carbide Cutters with handpiece shank for a broad field of application in oral surgery. The cutters have a particularly easy cutting and efficient saw tothing and a specifically designed twist angle. Further, the narrow working parts allow for an optimal view on the working area.

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Art.-No. ALV18

The Alveoplasty Surgical Kit contains the most important instruments for surgical contouring and modification of the alveolar ridge when unwanted unevenness of the jawbone occurs after tooth extraction or bone degeneration. In this case, the kit offers the best way to smooth the jawbone. The Alveoplasty Kit contains two initial burs and two Lindemann burs, which allow for the performance of precise, initial osteotomies, even in hard cortical bone. To perform the alveoplasty itself, five different carbide cutters and a round diamond are included. The round working parts optimally protect the soft tissue against injuries. Thus, the alveolar ridge is prepared carefully and precisely for implant insertion or denture restoration.

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The Implant Preparation Kit provides the easy and safe preparation of the implant site. The stop sleeves allow for a defined depth stop during the initial and pilot drilling. For the application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the pilot drill of 0.8 mm must be considered.

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AT A GLANCE

- Simple and safe preparation of the implant site
- Defined depth stop during initial and extension drilling
- Colored marking of the stop sleeves for easy differentiation
- Additional depth-marked instruments for optional use without stop sleeves

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
The Drill-Stop-Control system provides an easy and safe preparation of the implant site. The increasing diameters of the twist drills are perfectly coordinated with each other. The stop sleeves provide the defined depth stop during the pilot drilling and during the expansion of the implant site. For the application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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

The Drill-Stop-Control BLT provides an easy and safe preparation of the implant site for Bone Level Tapered implants of the Straumann® or Oktagon® implant system. The increasing diameters of the twist drills are perfectly coordinated with each other. The stop sleeves provide the defined depth stop during the pilot drilling and during the expansion of the implant site. For the application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.4 mm must be considered.

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
With the specially for myplant two developed Drill-Stop-Control system, the preparation of the implant site is particularly easy and safe, since the diameters of the twist drills are perfectly matched to the myplant two implant system. The stop sleeves are used both during the pilot and the extension drilling for the defined depth stop. The lengths of the stop sleeves are ideal for subcrestal placement of the implant. For application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra lengths of the instruments must be considered.

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

Diameters fit perfectly with the myplant two implant system.

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
The Guided-Drill-Stop system is used for the easy and safe preparation of the implant site. The special feature of the system are the specially developed twist drills with guiding tip, with which an intermediate step is added prior to the actual extension drilling. With this step, the upper part of the implant site is expanded to the desired diameter, while the lower part is used to guide the drill. This avoids deviation from the desired axis. In addition, the stop sleeves provide the defined depth stop during all drillings.


CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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

3D-Navigation-Control is a guiding sleeve system that is suitable for all commonly used implant systems and planning programs. Due to the included closed standard sleeves and the open guiding sleeves for the molar region, the system permits the use of surgical guide templates in the anterior region and even in tight molar situations. In addition, the system provides special tool requirements of flapless surgery. The 3D-Navigation-Control system contains a separate 3D-Laboratory Kit so the dental technician is able to install the sleeves in the template precisely.

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Art.-No. BNA00

© Dr. Hans-Joachim Nickenig

Inclus | content : 3D-Navigation-Control & 3D-Navigation-Control Professional

Ref.	HN010	HN011	HN012	HN013	HN014	HN015	HN016	HN017	HN018
Tige ¹	205	206	206	206	206	206	206	206	206
Taille ²	020	020	025	030	035	040	020	025	030
Longueur mm	16.0	19.0	19.0	19.0	19.0	19.0	22.0	22.0	22.0
* ³	2.0	2.0	2.5	3.0	3.5	4.0	2.0	2.5	3.0

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3D-Navigation-Control Professional

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
3D-Navigation-Control Professional is a guiding sleeve system that is suitable for all commonly used implant systems and planning programs. Due to the included closed standard sleeves and the open guiding sleeves for the molar region, the system permits the use of surgical guide templates in the anterior region and even in tight molar situations. In addition, the system provides special tool requirements of flapless surgery. The 3D-Navigation-Control Professional system contains a separate 3D-Laboratory Kit so the dental technician is able to install the sleeves in the template precisely.

The Professional system offers an additional, separate 3D-Fixation Kit, which allows for the (optional) fixation of the surgical guide template.

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

The Periimplantitis Kit contains eight perfectly matched tungsten carbide finishing burs in egg and flame shapes for the idea peri-implantitis treatment. Both shapes come in two sizes and each in standard and ultra-fine toothing. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access, deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counter-clockwise so the implant surface is smoothed.

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AT A GLANCE

- Perfectly matched tungsten carbide finishing burs for the effective, intraoral treatment of titanium
- XXL FG shanks facilitate the treatment of difficult to access implants
- Different shapes and sizes for every individual patient case
- Clean and smooth titanium surfaces

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
The Degranulation Kit contains four extra coarse, round diamond instruments with a diameters of 1.0, 2.5, 3.0 and 3.5 mm. The instruments have been specially developed to enable quick and easy bone debridement. With the help of these diamonds, degranulation tissue adhering to the bone can be removed quickly and easily without damaging the surrounding bone, so that subsequently inserted grafting material comes into direct contact with healthy bone. This makes the Degranulation Kit a crucial addition to the degranulation process.

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AT A GLANCE


- Extra coarse, round diamonds for quick and easy debridement
- Reliable removal of degranulation tissue without damaging the surrounding bone
- Four different diameters for use in all situations
- Extra long instrument necks for an optimal view of the operating area
- Optimal addition to the degranulation process

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MEISINGER's dPTFE membranes are made of high-density PTFE with a pore size of less than 0.3 μm . Thus, the non-resorbable membrane is insensitive to bacteria and allows a maximum of predictability and aesthetics. When using the dPTFE membranes, no primary wound closure is necessary. The exposed membrane can be removed with tweezers after at least 21 days without surgical intervention, no additional anesthesia is required. The hexagonal depressions are designed to increase the surface area, which leads to an increased stability of the membrane and favors cell adhesion without increasing the porosity. For best results, the side with the depressions should be placed in the direction of the soft tissue. The Meisinger dPTFE membranes can be cut to size as needed.

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AT A GLANCE

- Non-resorbable membrane made of high-density PTFE
- Maximum of predictability and aesthetics
- No primary wound closure necessary
- Non-surgical removal of the membrane
- Hexagonal depressions of the surface for high stability
- No disruption of wound healing due to early resorption

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MEISINGER's Titanium-Reinforced membranes are made of high density PTFE. In combination with titanium reinforcement, the membranes offer exceptional dimensional stability and optimum flexibility at the same time. Therefore, the membranes are particularly suitable for the treatment of 3- to 4-wall defects in the course of the alveolar ridge augmentation. The wide variety of membrane shapes and sizes covers a variety of different defects and minimizes the need for trimming. If necessary, however, the membranes can be optimally trimmed to the required shape. The titanium-reinforced membranes are available in thicknesses of 150 μm and 250 μm .

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AT A GLANCE


- Non-resorbable membrane made of high-density PTFE
- Fine Titanium Reinforcement from Titanium Grade I
- Optimum combination of flexibility and dimensional stability
- Treatment of 3- to 4-wall defects
- Wide range of different membrane shapes
- Availability of two different membrane thicknesses

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
Designed for narrow single-tooth extraction sites,
especially where one or more bony walls are missing.

Art.-No. MTi2501230-1 Art.-No. MTi2501230-2
Art.-No. MTi1501230-1 Art.-No. MTi1501230-2
12 mm x 30 mm 12 mm x 30 mm
(1 membrane / boîte) (2 membranes / boîte)

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
Designed for grafting posterior extraction sites and limited ridge
augmentation. Designed for large buccal defects.

Art.-No. MTi2502530-1 Art.-No. MTi2502530-2
Art.-No. MTi1502530-1 Art.-No. MTi1502530-2
25 mm x 30 mm 25 mm x 30 mm
(1 membrane / boîte) (2 membranes / boîte)

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Designed for grafting large bony defects, including ridge
augmentation.

Art.-No. MTi2503040-1 Art.-No. MTi2503040-2
Art.-No. MTi1503040-1 Art.-No. MTi1503040-2
30 mm x 40 mm 30 mm x 40 mm
(1 membrane / boîte) (2 membranes / boîte)

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Designed for the largest bony defects, including ridge
augmentation.

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Designed for single-tooth extraction sites, especially
where one or more bony walls are missing.

Art.-No. MTi2502438-1 Art.-No. MTi2502438-2
Art.-No. MTi1502438-1 Art.-No. MTi1502438-2
24 mm x 38 mm 24 mm x 38 mm
(1 membrane / boîte) (2 membranes / boîte)

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
Designed for bony defects between adjacent teeth,
including ridge augmentation.

Art.-No. MTi25038TC-1 Art.-No. MTi25038TC-2
Art.-No. MTi15038TC-1 Art.-No. MTi15038TC-2
38 mm x 38 mm 38 mm x 38 mm
(1 membrane / boîte) (2 membranes / boîte)

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Designed for large bony defects between adjacent
teeth, including ridge augmentation.


Art.-No. MTi25038D-1 Art.-No. MTi25038D-2
Art.-No. MTi15038D-1 Art.-No. MTi15038D-2
38 mm x 38 mm 38 mm x 38 mm
(1 membrane / boîte) (2 membranes / boîte)

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Designed for large bony defects, including distal exten-
sion of the posterior ridge.

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Designed to fit periodontal defects in the anterior.

Art.-No. MTi2501318-1 Art.-No. MTi2501318-2
13 mm x 18 mm 13 mm x 18 mm
(1 membrane / boîte) (2 membranes / boîte)

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Designed to fit periodontal defects in the posterior.

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
MEISINGER's porcine collagen membrane is a natural membrane made from porcine pericardium. The special manufacturing process ensures that the microporous, 3-layered architecture is preserved and no cross-linking chemicals are needed. This extraordinary structure leads to a very flexible membrane with particularly high tensile strength.

The porcine collagen membrane is characterized by its smooth handling, high flexibility and ideal adaptability. Once hydrated the membrane loses its original shape and becomes flexible, so that it optimally adapts to all defect surfaces. As the surface of the membrane is identical on both sides, both sides can be placed against the defect.

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AT A GLANCE

- Made of porcine pericardium
- Special, 3-layered architecture for high tensile strength
- Smooth handling
- Especially high flexibility and ideal adaption to the defect
- Both sides of the membrane can be placed against the defect
- Resorption time of about 26 weeks

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
MEISINGER's bovine collagen membrane is made from highly purified type 1 collagen derived from bovine achilles tendon. A great advantage of the bovine membranes is their ideal handling. They combine high flexibility with optimum stability. This makes them flexible, adaptable to any shape and easy to place. At the same time they are stable enough to hold their shape and not to collapse into the defect. The unique fiber orientation ensures maximum tensile strength. Thus, the membranes can be safely fixed by sutures or pins without tearing. Due to the long and predictable resorption time of 26 to 38 weeks, there is no risk of particle loss due to early resorption.

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AT A GLANCE

- Long resorption time of 26 to 38 weeks
- Ideal handling
- High flexibility for easy ridge preservation
- Optimum dimensional stability
- High tensile strength for safe fixing

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
The high quality of the MEISINGER PTFE suture allows for atraumatic wound closure after dental procedures. The high-quality needles are made of stainless steel and have a very high bending strength with a minimum of wear. Sharp needles are the result, even after repeated piercing. The constant cutting performance allows for a particularly tissue-saving work. The suture is made of 100% medical grade bioinert PTFE. Thanks to its smooth surface, the monofilament prevents the adhesion of bacteria to the thread and ensures a resistance-free gliding through the tissue. The PTFE is a soft, non-resorbable material with consistent tensile strength, allowing the user to work very precisely, planning with a defined healing time, and patient-friendly care.

Remarque: 12 Sutures / boîte, Longueur suture : 45 cm
 Note: 12 Sutures / box, Suture length: 45 cm

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The autoclavable Tecapro™ Storage Tray is cleverly designed to store not only fixation instruments, but also a variety of membrane fixation, bone fixation, and tenting screws needed in dental bone grafting.

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Up to 100 fixation screws – including bone fixation, membrane fixation, and tenting – may be stored in a built-in screw organizer dial designed and labeled for easy identification, easy storage and simple re-ordering.

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The stainless steel driver handle and cruciform driver blade fit together to make up a drive system that ensures easy pick-up of screws and stable transfer to the surgical site.

Composants individuels, Individual Components

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Blades work universally with all Pro-Fix™ screws.

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Pro-Fix™ Membrane Fixation Screws are designed as an attractive alternative to using tacks for membrane stabilization. Easy pick-up, solid stability of the screw during transfer to the surgical site, and easy placement make membrane fixation fast and easy.

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Pro-Fix™ Tenting Screws are designed with a self-drilling tip, polished neck and broader head to maintain space under resorbable and non-resorbable membranes in horizontal and vertical bone regeneration procedures.

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Pro-Fix™ Bone Fixation Screws are designed with finer pitched, self-tapping threads that give the screws greater clamping force while using less driver torque. The screws' threads are equipped with a cutting flute that allows for easier insertion into harder bone. The screws are placed into a 1.2 mm pre-drilled pilot hole.

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Further self-tapping Bone Fixation Screws

N° art.	Contenu	
packed à	Content	
PFB8	1.5 mm x 8.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 8.0 mm self-tapping Bone Fixation Screw	
PFB8-5	1.5 mm x 8.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 8.0 mm self-tapping Bone Fixation Screws	
PFB10	1.5 mm x 10.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 10.0 mm self-tapping Bone Fixation Screw	
PFB10-5	1.5 mm x 10.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 10.0 mm self-tapping Bone Fixation Screws	
PFB12	1.5 mm x 12.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 12.0 mm self-tapping Bone Fixation Screw	
PFB12-5	1.5 mm x 12.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 12.0 mm self-tapping Bone Fixation Screws	
PFB14	1.5 mm x 14.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 14.0 mm self-tapping Bone Fixation Screw	
PFB14-5	1.5 mm x 14.0 mm vis auto-taraudantes de fixation de greffons 1.5 mm x 14.0 mm self-tapping Bone Fixation Screws	

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NanoBone® Catalogue
separately available

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AT A GLANCE

Synthetic bone grafting material

NanoBone® contains absolutely no human or animal ingredients so no material-related contamination risks and alleviates ethical concerns for patient consent. High-tech production processes ensure consistent and excellent quality.

Complete remodelling

Thanks to its special structure, NanoBone® can be used for natural remodelling. Osteoclasts resorb the material as osteoblasts form new bone. Material resorption and bone formation are linked. This means predictable results and a reliable bone basis.

Controlled osteoinduction

NanoBone® is osteoconductive and osteoinductive. In other words, it is not simply a scaffold: it actively promotes bone formation.

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AT A GLANCE

The applicator


- "Ready to use" direct application from the applicator
- Easy, intuitive use
- Precise placement of the NanoBone® QD material
- Controlled and exact portioning
- Allows for small approaches to the application site (outer diameter of the applicator just 8 mm)
- Allows for retrograde filling of bone defects
- Difficult to access bone defects treated due to the long applicator

The material

- Shapeable and adaptable to the defect
- Strong cohesion of the material
- Sticks to the bone (adhesion)
- Stable also in aqueous environment and profusely bleeding wounds
- Controlled osteoinduction

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
Le fermoir avant s'enlève en se tournant. The front safety lock is released by turning and then removed.

NanoBone® QD (39% Silice / 61% Hydroxylapatite)		206 m ² surface interne/internal surface area	
N° art.	Contenu/Content	Utilisations recommandées/Recommended applications ¹	
23P074	3 x 0,25 ml	<ul style="list-style-type: none"> Extraction alvéolaire (dents antérieures) Augmentation autour d'implants Péri-implantites Defauts parodontaux Apicoectomie 	<ul style="list-style-type: none"> Extraction alveoli (anterior tooth) Augmentation around implants Peri-implantitis Periodontal defects Apicoectomy
23P075	3 x 0,5 ml	<ul style="list-style-type: none"> Extraction alvéolaire (dents antérieures) Augmentation autour d'implants Péri-implantites Defauts parodontaux Remplissage de kystes Apicoectomie 	<ul style="list-style-type: none"> Extraction alveoli Augmentation around implants Peri-implantitis Periodontal defects Filling cysts Apicoectomy
200070	1 x 1,0 ml	<ul style="list-style-type: none"> Extraction alvéolaire Apicoectomie Élévation du plancher sinusal (environ 0,6ml par implant) 	<ul style="list-style-type: none"> Extraction alveoli Apicoectomy Sinus floor elevation (ca. 0.6 ml per implant)
200071	1 x 2,5 ml	<ul style="list-style-type: none"> Élévation du plancher sinusal (environ 0,6ml) Remplissage de kystes 	<ul style="list-style-type: none"> Sinus floor elevation (ca. 0.6 ml per implant) Filling cysts

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
The MEISINGER Easy-Clean Scissors made of high-quality stainless steel can be easily and completely disassembled and put back together again. Residues below the hinge can be completely removed – for a residue-free processing*. For dismantling, the scissors only need to be opened at a right angle. Then the upper part of the scissors is lifted briefly. In addition, the scissors have a matte coating that is applied using PVD procedure (vacuum-based coating). Thanks to this black coating the products reach high hardness and durability. Furthermore, disturbing light reflections are reduced, and glare-free work is possible.

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AT A GLANCE

- Surgical scissors made of high-quality stainless steel
- Easy and complete dismantling for residue-free processing*
- Matte coating for high hardness and durability as well as glare-free work
- Scissors in two different versions

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
Instruments manuels / Manual instruments

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
This ratchet was designed to transfer turning moments up to 70 Ncm. May be disassembled for simple and safe cleaning.

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
Torque ratchet for working with defined torque. In blocking position, torques for both the insertion and the loosening of screws can be transmitted, e.g. for the insertion of implants.

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
Ratchet for working in torque function and in blocked function. In torque function the application takes place using the bending rod, from which torques between 10 and 70 Ncm can be read. In blocked function, the ratchet enables a greater torque to be transmitted, for example when inserting implants or loosening connections.

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MEISINGER offers one of the most extensive programs in the area of the Surgical Round Drills, Surgical Cutters, and Trephines. The optimized special construction of the instruments guarantee a safe application and fast achievement for the desired working result.

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Surgical Cutters of Tungsten Carbide

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The special, easy-cutting and efficient saw toothing and the specifically selected angle of twist allow a wide field of application of the instruments within the field of oral surgery.

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







Surgical round drill of tungsten carbide

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

165RF Acier inoxydable/ stainless steel

Ref.	Tige / Shank	Nr de ref / Ref.-No.	2	2
			L mm	
165RF	HP	330 104 408 297	023	7,0
	RA L	330 205 408 297	023	

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
Fraise boule en acier

Surgical round drill of stainless Steel

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


*Inclus dans le kit Saw Basic Kit 7150 * contained in the Saw Basic Kit 7150

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The MEISINGER gingiva trimmers are intended for use in dental mucosal surgery. They have been specially developed for the versatile and gentle treatment of gingival tissue. The working part consists of particularly high-quality and stable zirconia. The main application focus of the MEISINGER gingiva trimmers is the safe and gentle removal of unwanted granulation tissue and hyperplastic gingiva (papillectomy) - even in hard-to-reach areas. Furthermore, the trimmers are used in the enlargement of the sulcus after crown preparation, for the best possible exposition of the preparation margin and the creation of an exact impression. The gingiva trimmers are also particularly suitable for the tissue-saving exposure of covered healed implants and for use in periodontal treatment. The application takes place without cooling, so that the resulting rotational energy can be used for tissue modulation. Thanks to heat coagulation, bleeding tendency is reduced.

Instructions for use:

- Application without cooling
- Optimal speed: 300.000 - 450.000 rpm

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Stainless Steel Bur Blocks, non-corrosive, sterilizable

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Hager & Meisinger GmbH | Hansemannstr. 10 | 41468 Neuss | Germany