MEMBRANES & SUTURES

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Master Line by Dr. Istvan Urban Master Line par le Dr. Istvan Urban

Master-Core	12
Master-Mill	13
Master-Pin-Control	14-15

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Trephine Ejection Kit	16-17
Micro Screw System	18-19

Surgical Instruments Instruments chirurgicaux

Saw Basic Kit	20
Microsurgical Instruments Instruments de microchirurgie	21
Manual Instruments Instruments manuels	21

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
MEISINGER's porcine collagen membrane is a natural membrane made from porcine pericardium. The special manufacturing process ensures that the microporous, 3-layered architecture is preserved and no cross-linking chemicals are needed. This extraordinary structure leads to a very flexible membrane with particularly high tensile strength.

The porcine collagen membrane is characterized by its smooth handling, high flexibility and ideal adaptability. Once hydrated the membrane loses its original shape and becomes flexible, so that it optimally adapts to all defect surfaces. As the surface of the membrane is identical on both sides, both sides can be placed against the defect.

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Art.-No. MPC3040
30 mm x 40 mm
(1 membrane / box)

CONVINCING ADVANTAGES

- Made of porcine pericardium
- Special, 3-layered architecture for high tensile strength
- Smooth handling
- Especially high flexibility and ideal adaption to the defect
- Both sides of the membrane can be placed against the defect
- Resorption time of about 26 weeks

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
Membrane en collagène bovin

Type 1 bovine Collagen Membrane
Membrane en collagène bovin de type 1

MEISINGER's bovine collagen membrane is made from highly purified type 1 collagen derived from bovine achilles tendon. A great advantage of the bovine membranes is their ideal handling. They combine high flexibility with optimum stability. This makes them flexible, adaptable to any shape and easy to place. At the same time they are stable enough to hold their shape and not to collapse into the defect. The unique fiber orientation ensures maximum tensile strength. Thus, the membranes can be safely fixed by sutures or pins without tearing. Due to the long and predictable resorption time of 26 to 38 weeks, there is no risk of particle loss due to early resorption.

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Art.-No. MBC3040
30 mm x 40 mm
(2 membranes / box)

CONVINCING ADVANTAGES

- Long resorption time of 26 to 38 weeks
- Ideal handling
- High flexibility for easy ridge preservation
- Optimum dimensional stability
- High tensile strength for safe fixing

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PTFE: Polytetrafluoroethylene (PTFE) is comprised of a carbon chain with two fluorine atoms for every carbon atom. The complete fluorination of the carbon chain, along with the strength of the carbon-to-fluorine bonds, makes PTFE highly stable. This stability results in a synthetic polymer that is non-resorbable, biologically inert and chemically non-reactive, and therefore an ideal material for many medical device applications. In addition to its long history in the field of guided tissue regeneration (GTR), PTFE has been used for over 30 years in cardiovascular applications such as suture, vascular grafts and heart valves.

Expanded PTFE: PTFE as a biomaterial differs in porosity based on the amount of expansion applied during manufacturing. Heating PTFE and then applying force expands the material's microstructure to make expanded PTFE (ePTFE). Under scanning electron microscopy, we see a network of dense nodes connected by fibrils. As the nodes and fibrils are expanded, the porosity of the material continues to increase.

Expanded PTFE has a long history of success in GTR procedures, particularly in periodontics. However, the highly porous structure of ePTFE allows ingrowth of bacteria when the membrane is exposed in the mouth. Exposure results in high rates of infection and frequently requires early removal of the device. In addition, the highly porous structure allows soft tissue ingrowth, which complicates removal, often requiring sharp dissection and extensive surgery. Expanded PTFE must be completely buried and primary closure must be maintained to ensure predictability. While expanded PTFE is useful and quite predictable in deep, buried sites for guided tissue regeneration, there is currently no role for this material in extraction site grafting where exposure is likely.

Dense PTFE: Dense PTFE, also known as high-density PTFE or dPTFE, is manufactured to eliminate expansion of the nodes and fibrils, resulting in a micro-porous material that is impervious to bacteria while still allowing diffusion of gases and small

molecules. Dense PTFE was designed to withstand exposure in the oral environment, which represents an improvement to earlier versions of ePTFE in many applications, especially socket preservation where deliberate membrane exposure offers several advantages.

Upon implantation, dense PTFE is immediately coated with plasma proteins, facilitating cellular adhesion to the smooth, biocompatible surface. This cellular adhesion is observed to form a hermetic seal, providing resistance to migration of bacteria and epithelial cells around and under the membrane when it is exposed in the mouth. Plasma protein adsorption also facilitates diffusion of soluble organic molecules across the membrane. Removal of dense PTFE is simplified due to the lack of tissue ingrowth into the surface structure.

A textured dense PTFE is available. Texturing the membrane results in an increase in surface area and may increase the pull-out strength of the material through three dimensional attachment of soft tissue. The increased stability in the wound may result in less flap retraction and reduce risk of membrane movement and loosening. The primary advantage of dense PTFE is the ability to remain exposed in the mouth while protecting the underlying defect and bone graft. The membrane is soft, flexible and easy to handle. Primary closure is not required, and the membrane may be removed without additional surgery if exposed. If primary closure technique is used, the membrane may be easily removed through a small incision in a flapless technique.

Dense PTFE is also available with titanium reinforcement, which increases the stiffness of the material for use in defects where spacemaking is required. The embedded titanium framework allows the membrane to be shaped to fit a variety of defects without rebounding and provides additional stability in large, non-space-making osseous defects.

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


1. Preoperative view: To minimize trauma to the alveolar bone and to optimize the result of the ridge preservation, a minimally invasive atraumatic extraction technique should be used. It is recommended to use periostomes or surgical sectioning to minimize mechanical trauma to the thin cortical bone.
2. All soft tissue remnants should be gently removed with a sharp curettage. Especially with endodontically treated teeth special care should be taken to remove residual soft tissues at the apical extent of the socket. Bleeding from the socket walls should be noted. If necessary, decortication of the socket wall can be done with a round bur to improve blood supply. Bleeding points should be set to increase early vascularization and access to osteoprogenitor cells.
3. A subperiosteal pocket is created with a small periosteal elevator or curette, extending 3-5 mm beyond the socket margins (or defect margins) on the palatal and the facial aspect of the socket. In the esthetic zone, rather than incising and elevating the interdental papilla, it is left intact and undermined in a similar fashion. The dPTFE membrane will be tucked into this subperiosteal pocket.
4. Particulate augmentation material is placed into the socket with a syringe or curette. Ensure that the material is evenly distributed throughout the socket. In order to provide enough space for the ingrowth of blood vessels, the augmentation material should not be condensed too tightly.
5. The dPTFE membrane is trimmed to extend 3-5 mm beyond the socket walls and then tucked subperiosteally under the palatal flap, the facial flap and underneath the interdental papilla with a curette. The membrane should rest on bone 360° around the socket margins, if possible. Note that minimal flap reflection is necessary to stabilize the membrane.
6. Ensure that there are no folds or wrinkles in the membrane and that it lies

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passively over the socket. Remove any stray bone graft particles that may be present between the membrane and the flap. To prevent bacterial leakage under the membrane, take care to avoid puncturing the membrane, and do not overlap two adjacent membranes.

7. The membrane is further stabilized with a criss-cross PTFE suture. It is not recommended to suture through the membrane. Alternatively, interrupted sutures may be placed. The PTFE sutures, which cause minimal inflammatory response, are left in place for 10 to 14 days.
8. The membrane is removed, non-surgically, in 21 to 28 days. Sockets with missing walls may benefit from a longer time frame. Topical anesthetic is applied, then the membrane is grasped with a tissue forcep and removed with a gentle tug.
9. Studies have shown that by 21-28 days there is a dense, vascular connective tissue matrix in the socket and early osteogenesis is observed in the apical 2/3 of the socket.
10. Immediately following membrane removal, a dense, highly vascular, osteoid matrix is observed. The natural position of the gingival margin has been left intact because primary closure was not necessary. The dense PTFE membrane has contained the graft material and prevented epithelial migration into the socket. Adjacent gingival epithelium migrates across the osteoid matrix upon removal of the membrane.
11. The socket after 6 weeks. Keratinized gingiva is beginning to form over the grafted socket. The natural soft tissue architecture is preserved, including the interdental papillae. New bone is beginning to form in the socket. Over the next 6 to 10 weeks, increasing thickness of trabeculae and mineralization will result in load bearing bone suitable for implant placement.

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
MEISINGER's dPTFE membranes are made of high-density PTFE with a pore size of less than 0.3 μm . Thus, the non-absorbable membrane is insensitive to bacteria and allows a maximum of predictability and aesthetics. When using the dPTFE membranes, no primary wound closure is necessary. The exposed membrane can be removed with tweezers after at least 21 days without surgical intervention, no additional anesthesia is required. The hexagonal depressions are designed to increase the surface area, which leads to an increased stability of the membrane and favors cell adhesion without increasing the porosity. For best results, the side with the depressions should be placed in the direction of the soft tissue. The MEISINGER dPTFE membranes can be cut to size as needed.

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Art.-No. MDP2530
25 mm x 30 mm
(4 membranes / box)

CONVINCING ADVANTAGES

- Non-absorbable membrane made of high-density PTFE
- Maximum of predictability and aesthetics
- No primary wound closure necessary
- Non-surgical removal of the membrane
- Hexagonal depressions of the surface for high stability
- No disruption of wound healing due to early resorption

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The MEISINGER Titanium-Reinforced membranes are made of high density PTFE. In combination with titanium reinforcement, the membrane offers exceptional dimensional stability and optimum flexibility at the same time. Therefore, the membrane is particularly suitable for the treatment of 3- to 4-wall defects in the course of the alveolar ridge augmentation. The wide variety of membrane shapes and sizes covers a variety of different defects and minimizes the need for trimming. If necessary, however, the membranes can be optimally trimmed to the required shape. The titanium-reinforced membranes are available in thicknesses of 150 μm and 250 μm . The 150 micron membranes are 40% thinner than the 250 micron membranes.

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© Hager & Meisinger GmbH, Germany

CONVINCING ADVANTAGES

- Non-absorbable membrane made of high-density PTFE
- Fine Titanium Reinforcement from Titanium Grade I
- Optimum combination of flexibility and dimensional stability
- Treatment of 3- to 4-wall defects
- Wide range of different membrane shapes
- Availability of two different membrane thicknesses

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
MEISINGER PTFE Suture Suture PTFE

Soft monofilament Suture
Suture monofilament souple

The high quality of the MEISINGER PTFE suture allows for atraumatic wound closure after dental procedures. The high-quality needles are made of stainless steel and have a very high bending strength with a minimum of wear. Sharp needles are the result, even after repeated piercing. The constant cutting performance allows for a particularly tissue-saving work. The suture is made of 100% medical grade bioinert PTFE. Thanks to its smooth surface, the monofilament prevents the adhesion of bacteria to the thread and ensures a resistance-free gliding through the tissue. The PTFE is a soft, non-resorbable material with consistent tensile strength, allowing the user to work very precisely, planning with a defined healing time, and patient-friendly care.

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CONVINCING ADVANTAGES

- 100% medical grade, bioinert high-quality-PTFE
- Monofilament does not absorb bacteria
- Soft thread for patient-friendly care
- No memory effect for excellent handling and safe knotting
- Non-resorbable material for reliable wound closure

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Art. No. MS4013TP
Suture: USP 4-0

Needle: 13 mm, 1/2 Circle, Round-Bodied
Aiguille : 13 mm, 1/2 de cercle, Ronde

Art. No. MS4013
Suture: USP 4-0

Needle: 13 mm, 3/8 Circle, Reverse Cutting
Aiguille : 13 mm, 3/8 de cercle, Tranchant inversé

Art. No. MS4016
Suture: USP 4-0

Needle: 16 mm, 3/8 Circle, Reverse Cutting
Aiguille : 16 mm, 3/8 de cercle, Tranchant inversé

Art. No. MS5013
Suture: USP 5-0

Needle: 13 mm, 3/8 Circle, Reverse Cutting
Aiguille : 13 mm, 3/8 de cercle, Tranchant inversé

Art. No. MS5016
Suture: USP 5-0

Needle: 16 mm, 3/8 Circle, Reverse Cutting
Aiguille : 16 mm, 3/8 de cercle, Tranchant inversé

Master-Core Basic & Master-Core Professional

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The Bone Management® Master-Core Basic and Professional Systems are used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the systems contain trephines with different lengths and diameters. Due to the different lengths, the trephines have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trephines' working parts enable additional depth control during the extraction. The black coating of the trephines allows for glare-free work and the optimal visibility of the depth markings.

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The Bone Management® Master-Mill is a surgical bone mill used for crushing autologous bone blocks and thus preserving small autologous bone chips. With the help of these bone chips, bone defects of all shapes can be optimally filled. The individual components of the Master-Mill are made of high quality, surgical stainless steel. All parts feature an optimal fit and a stable geometry, allowing the Master-Mill to be easy to handle with ease of assembly and intuitive use.

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Art.-No. BKM00

Fig. BKM00 BKM00

Take a Look at Urban Regeneration Institute For Courses and Information www.implant.hu

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The Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The included pins consist of titanium alloy Grade 5. Furthermore, they have a groove which increases the surface of the pins. In addition, due to the groove and a unique head-design, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.

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
¹ 204=RA ² Largest working part diameter in 1/10 mm

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Art.-No. BMPBA


Master-Pin-Control Basic
10x Pins Art.-No. BMPBA

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The Trephine Ejection Kit was developed together with Prof. Dr. Fouad Khoury. It contains several unique ideas which turn this kit into something special: Two-part, internally cooled trephines allow for a safe and minimal invasive extraction of local bone cylinders. The inner cooling helps to prevent a damage of the bone tissue due to overheating, while the removable trephine working parts allow for an especially easy sampling of the cut out bone cylinders. Either with the aid of the ejection needle or the ejection instruments, which perfectly match with the diameter of the trephine working parts, the bone cylinders can be pushed out of the trephines. Furthermore, the included predrilling trephines allow for a safe application of the trephines. Due to the two-part trephines, the variety of different diameters and the optimal matching of the instruments, the Trephine Ejection Kit offers highest efficiency, flexibility and precision to the user for the extraction of bone cylinders.

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© Prof. Dr. Fouad Khoury

Micro Screw System & Micro Screw System Basic

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The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1.0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter. The Micro Screw Systems can be used for vertical, horizontal and three-dimensional bone augmentation of different defect sizes.

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
Saw Basic Kit

Selection of various Diamond Coated Saws for Oral Surgery
 Selection de diverses Scies diamantées pour chirurgie orale

The Saw Basic Kit provides a selection of diamond saws in the three most important diameters. They are used for the application in bone spreading or splitting operations. Due to the various sizes, the set offers the clinician the possibility of having a proper selection of these saws during surgery.

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

Manual Instruments Instruments manuels

Ref.	TL7	TL8	TL19
Name	Membrane Positioning Positionnement de membrane	Condenser Condenseur	Applicator Applicateur
Length in mm	170.0	181.0	182.0

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Hager & Meisinger GmbH
Hansemannstrasse 10
41468 Neuss | Germany
Phone: +49 2131 2012-303
Fax: +49 2131 2012-222
E-Mail: info@meisinger.de
Internet: www.meisinger.de

Meisinger France sàrl
92 cours Lafayette (Mamaworks)
CS 53515
69489 Lyon Cedex 3
France
Téléphone : +33.4.26.78.74.83
E-mail : commande@meisinger.fr

**MANUFACTURER
FABRIQUANT**

**MEMBRANES & SUTURES
MEMBRANES & SUTURES**

 Osteogenics
4620 71st Street | Building 78-79
Lubbock, TX 79424 | USA

**MASTER-LINE, BONE TRANSFER & OSTEOSYNTHESIS,
SURGICAL INSTRUMENTS
MASTER-LINE, TRANSFERT OSSEUX & OSTÉOSYNTÈSE,
INSTRUMENTS CHIRURGICAUX**

 Hager & Meisinger GmbH
Hansemannstrasse 10
41468 Neuss | Germany